## Cross-Checks

REF CI 107

### Directions For Use

#### INDICATIONS FOR USE

The SteriTec Cross-Checks multi-variable (ISO 11140-1, Class 4) chemical indicator strips are designed for use in steam sterilizers operating at 132°C to 134°C (270°F to 273°F) for 3 minutes or longer. When used as directed, the SteriTec Cross-Checks indicator strips give visible indication that sterilizing conditions were met.

#### COLOR CHANGE

During steam sterilization, the check mark at the end of the strip changes from white to black.

#### CRITICAL PARAMETERS *(in a standard hospital steam sterilizer)*

Steam sterilization cycle functioning at 132° – 134°C for 3 minutes or longer.

#### STATED VALUES *(As determined in a steam sterilization resistometer.)*

134° C – 2.5 minutes

132° C – 3 minutes

#### INSTRUCTIONS FOR USE

1. Place a SteriTec Cross-Checks strip in each pack, peel-pouch or tray for sterilization. Place strip at a location in the package considered most difficult for steam to reach. For larger packs, place strips at several locations in pack.
2. Sterilize package as usual in a steam sterilizer.
3. Following sterilization cycle, verify that indicator check marks on the ends of the strip have turned from white to black.
4. If indicator is light brown, it suggests inadequate sterilization conditions and item(s) must be re-sterilized with new indicator strip(s).

#### SAFETY PRECAUTIONS

- *If there is any doubt about sterility of item, it must be considered NOT sterile.*
- *Do not use any strips whose indicator marks are not white in color prior to exposure to steam sterilization processing.*
- *Cross-Checks indicator strips must NOT be used in place of a biological indicator, but can be used in conjunction with a standard infection control monitoring protocol such as that described by the AAMI Standards and Recommended practices: Good Hospital Practice: Sterilization and Sterility Assurance.*

#### STORAGE

SteriTec Cross-Checks steam sterilization indicator strips should be stored between 50°-100°F (10°- 38°C) and at normal humidity (10 - 90% R.H.) in the re-sealable bag provided. No special storage conditions are necessary after exposure to steam sterilization conditions.

#### EXPIRY DATE

The expiry date is printed on the product packaging.

#### LOT NUMBER

A unique identification code, **LOT**, is printed on each indicator strip and packaging labeling.

#### INTERFERING SUBSTANCES OR CONDITIONS

There are NO KNOWN INTERFERING SUBSTANCES OR CONDITIONS that could affect the intended use of the indicator or adversely affect the indicator performance.

#### RELEASE OF TOXIC SUBSTANCES

The indicator releases NO KNOWN TOXIC SUBSTANCES in sufficient quantities to cause either a health hazard, or a hazard to the intended properties of the product being sterilized before, during or after the sterilization process.

#### DECLARATION OF CONFORMITY

SteriTec Products Mfg. Co., Inc. declares that the *Cross-Checks* Indicator conforms to all applicable portions ISO 11140-1 Class 4 Multi-variable Indicators.

SteriTec Products Mfg. Co., Inc.  
74 Inverness Drive East  
Englewood, CO 80112  
(303) 660-4201 FAX (303) 660-4213

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**ID&CO S.r.l.**