Cross-Checks L

REF CI 111

Directions For Use

INDICATIONS FOR USE

SteriTec Cross-Checks L chemical multi-parameter (ISO 11140-1, Class 4) chemical indicator strips are designed for use in steam sterilizers operating at 132°C to 134°C (270°F to 273°F). When used as directed, the SteriTec Cross-Checks L indicator strips give visible indication that sterilizing conditions were met.

COLOR CHANGE

During steam sterilization, the check mark at the end of the strip changes from white to brown/black.

CRITICAL PARAMETERS

Steam sterilization cycle functioning at 132° – 134°C for 3 minutes or longer.

STATED VALUES (Determined in a steam sterilization resistometer)

134°C – 3 minutes

132°C – 3.5 minutes

INSTRUCTIONS FOR USE

1. Place a SteriTec Cross-Checks L strip in each pack, peel-pouch or tray for sterilization. Place strip at a location in the package considered most difficult for steam to reach. For larger packs, place strips at several locations in pack.
2. Sterilize package as usual in a steam sterilizer.
3. Following sterilization cycle, verify that indicator check mark on strip has turned from white to brown/black.
4. After exposure to steam sterilization, check mark must change color to a brown/black. If indicator is white or tan, it suggests inadequate sterilization conditions and item must be re-sterilized with new indicator strip(s). ***If there is any doubt about sterility of item, it must be considered NOT sterile.***

SAFETY PRECAUTIONS

- *Caution: Do not use any strips whose indicator marks are not white in color prior to exposure to steam sterilization processing.*
- *Caution: Cross-Checks indicator strips must NOT be used in place of a biological indicator, but can be used in conjunction with a standard infection control monitoring protocol such as that described by the AAMI Standards and Recommended practices: Good Hospital Practice: Sterilization and Sterility Assurance.*

STORAGE

SteriTec Cross-Checks L steam sterilization indicator strips should be stored between 50°-100°F (10°- 38°C) in the re-sealable bag provided. No special storage conditions are necessary after exposure to steam sterilization conditions.

EXPIRY DATE

The expiry date is printed on the product packaging.

LOT NUMBER

A unique identification code, [LOT], is printed on each indicator strip and packaging labeling.

INTERFERING SUBSTANCES OR CONDITIONS

There are NO KNOWN INTERFERING SUBSTANCES OR CONDITIONS that could affect the intended use of the indicator or adversely affect the indicator performance.

RELEASE OF TOXIC SUBSTANCES

The indicator releases NO KNOWN TOXIC SUBSTANCES in sufficient quantities to cause either a health hazard, or a hazard to the intended properties of the product being sterilized before, during or after the sterilization process.

DECLARATION OF CONFORMITY

SteriTec Products Mfg. Co., Inc. declares that the *Cross-Checks L* Indicator conforms to all applicable portions ISO 11140-1 Class 4 Multivariable Indicators.

SteriTec Products Mfg. Co., Inc. USA
74 Inverness Drive East
Englewood, CO 80112
(303) 660-4201 FAX (303) 660-4213

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ID&CO S.r.l.