

STERI-PAK KIT LF

Reusable Bowie-Dick Test Pack with blue to black ink technology

Catalog No. BD 109/S-T

- **IMPROVED INK CHEMISTRY AIDS INTERPRETATION**
Blue color remains if air leak is present.
- **HIGH STANDARDS**
- **CONTINUOUS INK DESIGN FROM CENTER TO EDGE**
Makes air leaks easier to see.
- **PERMANENT COLOR CHANGE**
Can be used as a permanent record when stored away from light sources.
- **ENVIRONMENTALLY SAFE**
Improved chemistry is "Lead-Free" and Toxic Heavy Metals free.



ORDERING INFORMATION

PART #	PACKAGED
BD 109/S-T	1 CASE/16 PACKS (Each pack 30 Uses)
BD 105	Holder Each (Sold separately)

SteriTec

PRODUCTS, INC.

74 INVERNESS DRIVE EAST
ENGLEWOOD, CO 80112 • USA
(303) 660-4201 • Fax: (303) 660-4213
Website: <http://www.steritecproducts.com>
E-Mail: info@steritecproducts.com