


STERI-PAK™ LF

BOWIE-DICK TEST SHEETS

Large size Bowie-Dick Test Sheets with blue to black ink technology

Catalog No. BD 110

- **IMPROVED INK CHEMISTRY AIDS INTERPRETATION**
Vivid blue color remains if air leak is present.
- **HIGH STANDARDS**
- **PERMANENT COLOR CHANGE**
Can be used as a permanent record when stored away from light source.
- **ENVIRONMENTALLY SAFE**
Improved chemistry is "Lead-Free" and Toxic Heavy Metals free.



ORDERING INFORMATION

<u>PART #</u>	<u>PACKAGING</u>
BD 110	6 BUNDLES OF 100/CASE

SteriTec

PRODUCTS, INC.

74 INVERNESS DRIVE EAST
ENGLEWOOD, CO 80112 • USA

(303) 660-4201 • Fax: (303) 660-4213
Website: <http://www.steritecproducts.com>
E-Mail: info@steritecproducts.com