


INTERPRETATION GUIDE

Green Card T- BD 120

UNEXPOSED


GREEN CARD T



**AIR REMOVAL &
STEAM PENETRATION TEST**

**CORRECT COLOR
CHANGE**


GREEN CARD T



**AIR REMOVAL &
STEAM PENETRATION TEST**

INSUFFICIENT

GREEN CARD T



**AIR REMOVAL &
STEAM PENETRATION TEST**



74 Inverness Drive East - Englewood, Colorado 80112

Phone: 303-660-4201 - Fax: 303-660-4213 - E-mail: info@steritecproducts.com -

Website: www.steritecproducts.com

MK-083 Rev. A