

DISPOSABLE BIOLOGICAL TEST PACK—STEAM

With Instant Readout Integrator

Catalog No. BI 103

- **TEST PACK PARALLELS PERFORMANCE OF AAMI BIOLOGICAL TEST PACKS**

Easy open pull tab feature.


- **CONTAINS INTEGRATOR CARD WHICH PREDICTS THE BIOLOGICAL OUTCOME**

Integrator card inserted in the test pack will not complete a color change until after biologicals have been killed. This allows immediate release of the sterilizer load while the biological verification is incubating. Card provides a permanent record of the cycle.

- **CONTAINS SELF-CONTAINED BIOLOGICAL INDICATOR**

- **SMALL, COMPACT EASY TO USE**

The SteriTec Biological Test Pack for steam sterilization with Instant Readout Integrator is about the size of a deck of cards. Easy to store and easy to Use.



ORDERING INFORMATION

<u>PART #</u>	<u>PACKAGING</u>
BI 103	30 PACKS/CASE WITH 5 EXTRA CONTROLS

SteriTec

PRODUCTS, INC.

599 TOPEKA WAY • SUITE 400

CASTLE ROCK, CO 80109 • USA

(303) 660-4201 • Fax: (303) 660-4213

Website: <http://www.steritecproducts.com>

E-Mail: info@steritecproducts.com