

# CROSS-CHECKS-EO

## INTEGRATORS

An economical, inside-the-pack or pouch EO Integrator to check the combination of time, temperature, humidity and EO gas concentration.

Catalog No. CI 106

- SURPASSES PERFORMANCE OF BIOLOGICAL SPORE STRIPS
- ACCURATE

	TEST DATA	
	CROSS -CHECK EO	BIOLOGICAL SPORE STRIPS
TEMP<30 min.	FAIL	PASS
TEMP<110°F	FAIL	PASS
GAS<200 MG/L	FAIL	PASS
RH<20% RH	FAIL	PASS

Requires all parameters for EO gas sterilization be met. Distinctive yellow to brown color change is not reversible.

- **PLASTIC LAMINATION**  
Plastic lamination protects sterilized item from potential ink transfer
- **EASY TO INTERPRET**  
Color reference arrow printed on each indicator.

### ORDERING INFORMATION

PART #	PACKAGING
CI 106	24 BAGS OF 250/CASE



# SteriTec

### PRODUCTS, INC.

599 TOPEKA WAY • SUITE 400  
 CASTLE ROCK, CO 80109 • USA  
 (303) 660-4201 • Fax: (303) 660-4213  
 Website: <http://www.steritecproducts.com>  
 E-Mail: [info@steritecproducts.com](mailto:info@steritecproducts.com)