

CROSS-CHECKS

CHEMICAL INDICATORS

An economical, inside-the-pack indicator to check proper time and temperature in a pre-vacuum steam sterilizer operating at 132°C (270°F) U.S., 134°C (273°F) International

Catalog No. CI 134, CI 107,
CI 111

- **SURPASSES PERFORMANCE OF BIOLOGICAL SPORE STRIPS**
- **ACCURATE**
Requires all parameters for steam sterilization be met. Distinctive white to black or blue to black for CI 134. Color change is not reversible.
- **PLASTIC LAMINATION**
Short version (CI 111) Indicator is plastic laminated to protect sterilized items from potential ink transfer.
- **EASY TO INTERPRET**
Color reference arrow printed on each indicator.

ORDERING INFORMATION

<u>PART #</u>	<u>PACKAGING</u>
CI 107	16 BAGS OF 250/CASE
CI 111	24 BAGS OF 250/CASE
CI 134	24 BAGS OF 250/CASE



SteriTec

PRODUCTS, INC.

74 INVERNESS DRIVE EAST

ENGLEWOOD, CO 80112 • USA

(303) 660-4201 • Fax: (303) 660-4213

Website: <http://www.steritecproducts.com>

E-Mail: info@steritecproducts.com