

INTERPRETATION GUIDE

Cross Checks - CI 107

UNEXPOSED

	BROWN/MARRONE/CAFÉ = UNDEREXPOSED/INSUFFICIENTE/INSUFICIENTE BLACK/NERO/NEGRO = CORRECT/CORRETTO/CORRECTO	CROSS-CHECKS STEAM STERILIZATION MONITORS <i>SteriTec</i> LOT 090301	STEAM ISO 11140-1 CLASS/CLASSE/CLASE 4 MULTI-VARIABLE INDICATOR INDICATORE MULTIPARAMETRO/INDICADORA	
---	--	---	--	---

CORRECT COLOR CHANGE

	BROWN/MARRONE/CAFÉ = UNDEREXPOSED/INSUFFICIENTE/INSUFICIENTE BLACK/NERO/NEGRO = CORRECT/CORRETTO/CORRECTO	CROSS-CHECKS STEAM STERILIZATION MONITORS <i>SteriTec</i> LOT 090301	STEAM ISO 11140-1 CLASS/CLASSE/CLASE 4 MULTI-VARIABLE INDICATOR INDICATORE MULTIPARAMETRO/INDICADORA	
---	--	---	--	---

FAIL

	BROWN/MARRONE/CAFÉ = UNDEREXPOSED/INSUFFICIENTE/INSUFICIENTE BLACK/NERO/NEGRO = CORRECT/CORRETTO/CORRECTO	CROSS-CHECKS STEAM STERILIZATION MONITORS <i>SteriTec</i> LOT 090301	STEAM ISO 11140-1 CLASS/CLASSE/CLASE 4 MULTI-VARIABLE INDICATOR INDICATORE MULTIPARAMETRO/INDICADORA	
---	--	---	--	---



74 Inverness Drive East - Englewood, Colorado 80112

Phone: 303-660-4201 - Fax: 303-660-4213 - E-mail: info@steritecproducts.com -

Website: www.steritecproducts.com

MK-087 Rev. A