

STEAM CHARTS

CHEMICAL INTEGRATORS 121 °C / 134 °C

For use in testing steam sterilizers operating at 121 °C or at 134 °C

Catalog No. CI 108, CI 109,
CI 108/121-134

- **TEMPERATURE SPECIFIC**
Three different indicators.
CI 108, for testing at 121 °C.
CI 109, for testing at 134 °C
CI 108/121-134, for both 121 °C & 134 °C
- **PERFORMS EQUIVALENTLY TO BIOLOGICAL INDICATORS AT 121 °C AND AT 134 °C.**
- **MULTIPLE COLOR CHANGE INKS**
First spot tells if the item has been processed in a steam sterilizer.

Second spot of 134 °C indicator (third spot of 121 °C indicator) tells if correct conditions for sterilization have been met.

Center dot on all indicators shows excessive sterilization which can be detrimental to instruments and wastes energy.
- **COLOR CHANGE IS PERMANENT**
Color change from purple to green is not reversible.
- **PLASTIC LAMINATION**
Plastic lamination eliminates potential for ink transfer.

ORDERING INFORMATION

PART #	PACKAGING
CI 108	20 BAGS OF 250/CASE
CI 109	20 BAGS OF 250/CASE
CI 108/121-134	20 BAGS OF 250/CASE



CI 108



CI 109



CI 108/121-134

SteriTec

PRODUCTS, INC.

599 TOPEKA WAY • SUITE 400
CASTLE ROCK, CO 80109 • USA
(303) 660-4201 • Fax: (303) 660-4213
Website: <http://www.steritecproducts.com>
E-Mail: info@steritecproducts.com