

DRY HEAT INDICATOR LABELS

Catalog No. CI 113


- **DESIGNED FOR USE IN DRY HEAT OPERATING AT 160°C**
- **GIVES VISIBLE CONFIRMATION THAT THE ITEM HAS BEEN EXPOSED TO DRY HEAT**
- **THE GREEN INDICATOR BAR TURNS TO DARK BROWN/BLACK**

Complete color change occurs in 40 minutes.


ORDERING INFORMATION

<u>PART #</u>	<u>PACKAGING</u>
CI 113	10 BOXES OF 1000/CASE

UNEXPOSED



PASS



SteriTec

PRODUCTS, INC.

74 INVERNESS DRIVE EAST

ENGLEWOOD, CO 80112 • USA

(303) 660-4201 • Fax: (303) 660-4213

Website: <http://www.steritecproducts.com>

E-Mail: info@steritecproducts.com