



Product Technical Data Sheet

Form No.: FTP-009 Rev.: F

Product Name: SteriTec Dry Heat Indicator Label

Product Type: Process Indicator

Reorder Number: CI 113

Description: The SteriTec Dry Heat Indicator Labels are designed for use in a dry heat sterilizer operating at 160°C. The label provides a clear identification of processed versus unprocessed items. When used as directed, the Dry Heat Label gives visible confirmation that the item was exposed to a dry heat process.

Performance Characteristics:

When exposed to dry heat conditions, the indicator bar changes color.

Color Change: Green to Dark Brown/Black

Complete Color Change Time: 10 minutes +/- 2 minutes at 160°C.

Indicator Stability: Excellent before and after sterilization

Physical Characteristics:

Dimensions: 0.75" W X 1.14" L

Indicator Substrate Material: Removable adhesive label

Quantity / Box: 1000

Boxes / Case: 10

Selling Unit: Case

Environmental Limitations:

Storage: 50° - 100°F (10° - 38°C)

Shelf Life: The expiry date is printed on the product packaging.

FDA 510(k) Release to Market Date: April 2, 1997 - K953122

Standards Classifications:

ANSI/AAMI ST79, Class 4

ISO 11140-1, Class 4