

CROSS-CHECKS P

Chemical Indicator for Sterrad® Plasma Sterilizers

Catalog No. CI 115

- **AN ECONOMICAL PROCESS INDICATOR FOR MONITORING PLASMA STERILIZATION CYCLES.**
- **DISTINCT COLOR CHANGE**
From blue to pink.
- **STABLE COLOR**
Indicator ink will not fade in storage or under fluorescent lights. Completed color change is not reversible.
- **EASY TO INTERPRET**
Color reference arrows printed on each indicator. Completed color change of check marks should match or be lighter than arrows.



ORDERING INFORMATION

| PART # | PACKAGING |
|--------|---------------------|
| CI 115 | 24 BAGS OF 250/CASE |

SteriTec

PRODUCTS, INC.

74 INVERNESS DRIVE EAST
ENGLEWOOD, CO 80112 • USA

(303) 660-4201 • Fax: (303) 660-4213
Website: <http://www.steritecproducts.com>
E-Mail: info@steritecproducts.com