

CROSS-CHECKS F

Chemical Indicator for Steam-Formaldehyde sterilizers

Catalog No. CI 118

- **MONITORS STEAM-FORMALDEHYDE PROCESS**

An economical process indicator to monitor sterilization cycles in LTSF, Low Temperature Steam-Formaldehyde sterilizers.

- **DISTINCT COLOR CHANGE**


From purple to green

- **PERMANENT COLOR CHANGE**

Color change from purple to green is not reversible. Indicator will not fade in sunlight or change prematurely in storage. Can be used as permanent record.

- **HIGH STANDARDS**

Meets ISO 11140-1 standards for low temperature Steam-Formaldehyde process indicators.



ORDERING INFORMATION

PART #	PACKAGING
CI 118	24 BAGS OF 250/CASE

SteriTec

PRODUCTS, INC.

74 INVERNESS DRIVE EAST

ENGLEWOOD, CO 80112 • USA

(303) 660-4201 • Fax: (303) 660-4213

Website: <http://www.steritecproducts.com>

E-Mail: info@steritecproducts.com