



# Product Technical Data Sheet

Form No.: FTP-056 Rev.: F

**Product Name:** SteriTec Cross-Checks F

**Product Type:** Formaldehyde Indicator

**SteriTec Reorder Number:** CI 118

**Description:** The SteriTec Cross-Checks F was designed for use in Low Temperature Steam/Formaldehyde sterilizers. The indicator check mark will change from purple to green to give visual confirmation of exposure to sterilization conditions.

**Performance Characteristics:**

The indicator check mark will change color when exposed to steam/formaldehyde sterilization conditions.

**Color Change:** Purple to Green

**Indicator Stability:** Excellent before and after sterilization exposure

**Physical Characteristics:**

**Background Print Color:** Text = Black, Graphic = Green (Light 7489U)

**Dimensions:** 3.8125" X .625"

**Substrate Material:** Polystyrene

**Units / Bag:** 250

**Bags/Case:** 24

**Selling Unit:** Case

**Environmental Limitations:**

**Storage:** Room Temperature 10° – 38°C (50° - 100°F) and 10 - 90% relative humidity

**Shelf Life:** The expiry date is printed on the product packaging.

**SteriTec FDA 510(k) Release to Market Date:** Not Released. Not for USA distribution.

**Standards Classifications:**

ISO 11140-1