

# AUTOCLAVE TAPE

Steam/Chemical Vapor Sterilization Tape, Blue Steam Sterilization Tape, EO Sterilization Indicator Tape

Catalog No. CI 122, CI 124, CI 123

- **TAPES ARE STERILIZATION PACKAGE SEALING TAPES WITH A PROCESS INDICATOR PRINTED ON THE TAPE.**

Steam Tape indicator turns from white to dark brown/black when exposed to steam sterilization conditions.

EO Gas process indicator turns from brown to green when exposed to ethylene oxide sterilization conditions.

- **TAPE ADHERES TO ITEM THROUGHOUT THE STERILIZATION PROCESS.**

- **SHELF LIFE**

Steam - 3 years from manufacture date  
EO - 3 years from manufacture date

- **COMPLETE COLOR CHANGE**

Steam 121°C - 10 Minutes  
Steam 132°C - 134°C - 2 minutes  
Chemical Vapor 132°C - 20 Minutes  
EO Tape - 45 minutes +/- 5 minutes

- **60 YARDS (+/- 1 FT) PER ROLL**

- **STORE AT 23°C +/- 5°C, AWAY FROM LIGHT**



Steam Tape



Blue Tape



EO Tape

## ORDERING INFORMATION

PART #	PACKAGING
CI 122/.5	1/2" Steam Tape 72 rolls/CASE
CI 122/.75	3/4" Steam Tape 48 rolls/CASE
CI 122/1	1" Steam Tape 36 rolls/CASE
CI 124/.75	3/4" EO Tape 48 rolls/CASE
CI 124/1	1" EO Tape 36 rolls/CASE
CI 123/1	1" Blue Steam Tape 36 rolls/CASE

# SteriTec

## PRODUCTS, INC.

74 INVERNESS DRIVE EAST  
ENGLEWOOD, CO 80112 • USA

(303) 660-4201 • Fax: (303) 660-4213

Website: <http://www.steritecproducts.com>

E-Mail: [info@steritecproducts.com](mailto:info@steritecproducts.com)