



Product Technical Data Sheet

Form No.: FTP-167 Rev.: C

Product Name: SteriTec Dry Checks Dry Heat Indicator Strip

Product Type: Process Indicator

Reorder Number: CI 136

Description: The SteriTec Dry Checks Dry Heat Indicator Strips are designed for use in a dry heat sterilizer operating at 160°C. The indicator provides a clear identification of processed versus unprocessed items. When used as directed, the Dry Checks Indicator gives visible confirmation that the item was exposed to a dry heat process.

Performance Characteristics:

When exposed to dry heat conditions, the indicator bar changes color.

Color Change: Green to Black

Complete Color Change Time: 7 minutes +/- 2.5 minutes at 160°C

Indicator Stability: Excellent before and after sterilization

Physical Characteristics:

Dimensions: 0.75" W X 3.8" L

Indicator Substrate Material: 80 lb Cover

Quantity / Box: 250

Boxes/Case: 20

Selling Unit: Case

Environmental Limitations:

Storage: 50° - 100°F (10° - 38°C)

Shelf Life: The expiry date is printed on the product packaging.

FDA 510(k) Release to Market Date: K 952122, April 2, 1997

Standards Classifications:

ISO 11140-1, Class 4