



Product Technical Data Sheet

Form No.: FTP-170 Rev.: C

Product Name: SteriTec International Container Card

Product Type: Process Indicator

SteriTec Reorder Number: CI 138

Description: The SteriTec International Container Card is a steam sterilization process indicator designed for use with various types of sterilization containers. The Container Cards will function in sterilizers that operate at 121°C or 134°C. It provides a visible confirmation that the instrument container being sterilized has been exposed to steam sterilization processing.

Performance Characteristics:

When exposed to steam sterilization conditions, the bar changes color.

Color Change: White to Black

Complete Color Change Time: 2 minutes +/- 30 seconds at 121°C
3.5 minutes +/-30 seconds at 134°C

Indicator Stability: Excellent Before and After Exposure

Physical Characteristics:

The card can fit in various container card-holders by folding or not folding the top sizing flap to the appropriate size.

Background Print Color: Black/Gray

Dimensions: 1.531" W X 2.75" L or with flap folded down 1.4" W X 2.75" L

Substrate Material: 80 lb Cover

Units / Box: 500

Boxes/Case: 6

Selling Unit: Case

Environmental Limitations:

Storage: 50° – 100°F (10° – 38°C) at 30 - 60% humidity

Shelf Life: The expiry date is printed on the product packaging.

SteriTec FDA 510(k) Release to Market Date: None. Not released for sale in the United States.

Standards Classifications:

ANSI/AAMI ST79, Class 1

ISO 11140-1, Class 1