

CROSS-CHECKS LS

Catalog No. CI 139

- **CLASS 4 INDICATOR FOR ROUTINE PROCESSING.**

3 Minutes @ 134°C

12 Minutes @ 121°C

- **PLASTIC LAMINATION PROTECTS ITEMS FROM POTENTIAL INK TRANSFER.**
- **REQUIRES ALL PARAMETERS OF STEAM STERILIZATION ARE MET, NOT JUST HEAT.**
- **WILL NOT REVERSE COLOR AFTER EXPOSURE, VERY STABLE CHEMISTRY.**
- **TURNS BROWN/BLACK IN DRY (SUPERHEATED) STEAM.**
- **PACKAGED IN RE-SEALABLE FOIL POUCHES TO RETAIN FRESHNESS.**

ORDERING INFORMATION

<u>PART #</u>	<u>PACKAGING</u>
CI 139	20 BAGS OF 250/CASE

UNEXPOSED



PASS



FAIL



SteriTec

PRODUCTS, INC.

74 INVERNESS DRIVE EAST

ENGLEWOOD, CO 80112 • USA

(303) 660-4201 • Fax: (303) 660-4213

Website: <http://www.steritecproducts.com>

E-Mail: info@steritecproducts.com