



Product Technical Data Sheet

Form No.: FTP-198 Rev.: D

Product Name: SteriTec Wash-Checks Holder

Product Type: Washer Efficiency Monitor Accessory

SteriTec Reorder Number: WC 102 & WC 102/R

Description: The Wash-Checks Holder acts as a holding device for use with the Wash-Checks monitor. It allows for easy placement of the Wash-Checks Monitor throughout the washer chamber. The Wash-Checks Monitor (WC 101) is an easy to use single use prepared test for determining the effectiveness of the wash phase of a hospital washer/disinfector.

Performance Characteristics:

When the Wash-Checks Monitor is placed in the Holder the test soil is half obscured, thus simulating the joint of a surgical instrument. The removal of all test soil indicates effective cleaning.

Physical Characteristics:

The Wash-Checks Holder consists of a stainless steel device designed for easy placement on washer trays and shelves. The Holder has an outer section that partially covers the inserted Wash-Checks Monitor during use.

Dimensions: 3.50" L X 1.5" W X .6"

Substrate Material: Electro Polished Stainless Steel

Units / Box: 1

Boxes / Case: 6

Selling Unit: Case

Environmental Limitations: None

Shelf Life: Not Applicable

SteriTec FDA 510(k) Release to Market Date: Not Applicable

Standards Classification: N/A